

HL7

3.24.2.0 | User Guide



intelerad

[Intelerad.com](https://www.intelerad.com)

COPYRIGHT

© 2003-2024 DICOM Grid, Inc. dba Ambra Health. All Rights Reserved.

No portion of the contents of this publication may be reproduced or transmitted in any form or by any means without the express written permission of DICOM Grid, Inc. dba Ambra Health.

Information in this documentation is subject to change without notice and does not represent a commitment on the part of DICOM Grid, Inc. dba Ambra Health. The software described in this document is furnished under a license agreement or nondisclosure agreement and can be used and copied only in accordance with the terms of the agreement. No part of this document may be reproduced, or transmitted in any form or by any means, electronic or mechanical, including photocopying or recording, for any purpose, without the express written permission of DICOM Grid, Inc. dba Ambra Health.

Unless otherwise noted, names and clinical data used herein are fictitious. No association with any real person or clinical event is intended or should be inferred.

Except as expressly provided otherwise in writing, the information provided in this document is provided AS IS, without any condition or warranty, whether written, oral, express, implied, legal, or statutory. No condition or warranty of merchantability or fitness for a particular purpose applies to anything provided by Intelrad. Without limiting the foregoing, Intelrad nor any of its suppliers warrants this documentation is or will be accurate, error free or meets or will meet user requirements.

TRADEMARKS

Intelrad®, Intelrad Medical Systems, and the Intelrad logo, IntelePACS, InteleViewer™, Enterprise Viewer™, InteleOne®, InteleArchive™, InteleShare VNA™, InteleShare™, InteleShare Research™, InteleOrchestrator™, InteleHeart™, InteleScreen™, and InSight PACS™ are either registered trademarks or trademarks of Intelrad Medical Systems Incorporated or its affiliates.

THIRD-PARTY TRADEMARKS

Adobe, Acrobat, and Reader are either registered trademarks or trademarks of Adobe Systems Incorporated in the United States and/or other countries. Firefox is a registered trademark of Mozilla Foundation in the United States and other countries. Google Chrome browser is a registered trademark of Google Inc. in the United States and other countries. Internet Explorer is a registered trademark of Microsoft Corporation in the United States and other countries. Intel, Pentium, Pentium II Xeon, and Pentium III Xeon are trademarks or registered trademarks of Intel Corporation or its subsidiaries in the United States and other countries. Microsoft and Windows are either registered trademarks or trademarks of Microsoft Corporation in the United States and/or other countries. NVIDIA is a trademark or registered trademark of NVIDIA Corporation in the United States and/or other countries. Oracle and Java are registered trademarks of Oracle Corporation and or its affiliates. PowerScribe 360 is a trademark or registered trademark of Nuance Communications Inc. or its affiliates in the United States and/or other countries. Sun, Sun Microsystems, and Java are trademarks or registered trademarks of Sun Microsystems, Inc. in the United States and other countries. Safari, Mac, and OS X are trademarks of Apple Inc., registered in the United States and other countries.

All other brand names, product names, or trademarks belong to their respective holders.

INTENDED USE

InteleShare PACS software is intended for use as a primary diagnostic and analysis tool for diagnostic images for hospitals, imaging centers, radiologists, reading practices and any user who requires and is granted access to patient image, demographic and report information.

InteleShare ProViewer, a component of InteleShare PACS, displays, modifies and manages diagnostic quality DICOM images including 3D visualization and reordering functionality.

Lossy compressed mammographic images and digitized film screen images must not be reviewed for primary image diagnosis or image interpretations. Mammographic images may only be viewed using cleared monitors intended for mammography display.

Not intended for diagnostic use on mobile devices.

CONTRAINDICATIONS—None.

Caution: Federal law restricts this device to sale by or on the order of a physician. Fluency in English is required as a prerequisite for operating the Ambra PACS product.

This system does not replace the education, skill, and judgment of properly trained medical practitioners. Only properly trained and qualified individuals shall have access to and use InteleShare PACS and must know of its functionality, capabilities and limitations. Typical users of this system are trained health professionals, physicians, nurses, and technologists.

Downloaded Images, Workstations and Isolated Installs: You and your users must maintain InteleShare PACS with the most current versions, including available updates and upgrades. Delaying or refusing updates or upgrades following a recall may result in a non-compliant InteleSharePACS.

SAFETY ISSUES: InteleSharePACS is a medical device, and as such, must meet medical device safety and effectiveness requirements imposed by national regulations. Any unmonitored or unconnected use of InteleSharePACS, or use of InteleSharePACS without a valid right may put the health and safety of patients at risk as you will not be advised of the availability of any software patch, bug fix, update or upgrade nor will be informed of Field Safety Notices, Medical Device Recalls or Advisory Notices related to InteleSharePACS. Client and authorized users must consult national regulatory site(s) to be informed of Field Safety Notices, Medical Device Recalls or Advisory Notices related to InteleSharePACS. Intelrad does not have access to authorized users systems to implement corrections to prevent (or correct) occurrences of patient safety issues. You are responsible to flow down recall and patient safety information to your users. The user of the medical device must report any serious incident that has occurred in relation to the medical device to the manufacturer (Intelrad) and the competent authority having jurisdiction in their locale.

Referring Physicians Use: Images for authorized referring physicians may not be of diagnosis quality and should not be used for diagnostic purposes.



DICOM Grid, Inc. dba Ambra Health
305 Church at N Hills St, 6th Floor,
Raleigh, North Carolina 27609, USA



InteleShare UDI is: **+AMBRAPACS0/\$73.24.2.0D**

DECLARATION OF CONFORMITY

Declaration of Conformity is issued under the sole responsibility of the DICOM Grid, Inc. dba Ambra Health.

Title:	HL7 User Guide
Software version:	3.24.2.0 (Multiple Patient Multiple Use)
Date:	2024-05-14
Part number:	AMHL7EN3.24.2.0UG-O Issue 001

TABLE OF CONTENTS

InteleShare HL7 Guide	4
Introduction	5
Browser Requirements	6
Incoming HL7 Messages	7
ORU^R01 HL7 message	7
ORM^O01 HL7 message	8
Configuration	10
ADT^A08, ADT^A31, ADT^A04, ADT^A34 and SIU messages	12
Outgoing HL7 Messages	14
ORU^R01 HL7 messages	14
ORM^O01 HL7 messages	15
Message Segments and Fields	16
Recent Updates	17
Obtaining Printed Documentation	18
Contacting Intelrad Technical Support	19

3

INTELESHARE HL7 GUIDE

This guide provides an overview of IntelShare HL7 functionality. HL7 is one of the standard ways of passing medical information between systems.

4

INTRODUCTION

This guide provides an overview of the IntelShare HL7 functionality. HL7 is one of the standard ways of passing medical information between systems. IntelShare can both receive (incoming) and send (outgoing) HL7 messages. In terms of receiving, IntelShare supports orders (ORM^O01), reports (ORU^R01), and admit, transfer, discharge messages (ADT^A08 or ADT^A31). In terms of sending, IntelShare supports HL7 message types like orders (ORM^O01), reports (ORU^R01), and patient demographic queries (QRY^A19).

This guide will discuss the following topics:

- Incoming (to IntelShare) messages and workflows
- Outgoing (from IntelShare) messages and workflows
- Sample messages
- HL7 message segments and fields

IntelShare Gateway & Account Configuration

IntelShare transmits (sends and receives) HL7 messages via IntelShare Gateway. IntelShare Gateway is a piece of software that sits on a Microsoft® Windows® machine (VM or physical appliance) in the sending or target network. IntelShare also supports a cloud Gateway if preferred.

5

BROWSER REQUIREMENTS

Microsoft Internet Explorer 9 or later, Microsoft Edge, Apple Safari, Google Chrome, Firefox.

6

INCOMING HL7 MESSAGES

IntelSharecomplies with HL7 2.3. The messages IntelSharereceives as well as the key fields for these messages are outlined below.

ORU^R01 HL7 message

An ORU^R01 message is an unsolicited transmission of an observation. Colloquially, we refer to ORU^R01s as reports. It is highly useful to store the report with the imaging. IntelShare receives the ORU^R01 and ties the report to the study, typically using the accession number as the linking key. In addition, IntelShare uses HTML that IntelShare configures to make the ORU^R01 appear as a formatted report. These are the incoming ORU workflows that IntelShare supports:

Site sends an ORU^R01 to IntelShare

Upon receipt, IntelShare automatically queries and retrieves the study from the PACS or VNA. Upon receipt of the complete package (report and study), IntelShare automatically shares the package with the referring physician.

Site sends an ORU^R01 to IntelShare after the study already exists in IntelShare

As soon as the ORU^R01 is received, the report appears attached to the study. Those with the appropriate permissions can search for, find and review the study and report at any time.

Below is an example ORU^R01 message. The key fields are the accession number and the patient MRN. These are the linking keys between the report and the study. The

accession number is generally pulled from OBR1 or OBR2 and the MRN is generally pulled from PID2 or PID3 (see red text in sample message below). The field must contain only the accession number. It cannot contain any ^ or sub-fields. Although OBR2 or OBR3 is typically used, any segment and field in the ORU^R01 message can be defined as the field to link the report with the study.

See red text in sample message below.

```
MSH|^~\&|EMR|EMR|||20040111181023|1148|ORU^R01|487|D|2.3|| PID|1||1  
43560^^^^EPI||ECKLES^SAMUEL^A.^MR.^||19730211|M|ECKLES^SAMUEL^A^^^||4  
2 NORTH AVE^^CAMBR|InteleShareE^MA^02138^US^^^DN| DN|(857)721-3255|  
(857)721-3255|||1111246~11480172~7734 |114-80-2122|||^NEWYORK^^^
```

```
PD1|||FRANKIN NEAR EAST LOCATION (CENTRAL)^537| Outgoing Imaging Orders  
and Results Interface Technical Specification 10  
PV1|||^EMR^^^^^|||522320|||V ORC||437951^EPC|2223^EP  
C||Final||^200312071730^^^||200401211810|1152^KARIS^LUCAS^^^||757^SP  
ETZLER^CHRISTIAN^J^^^|467^^^537^^^|(857)259-9167||  
OBR|1|437951|2255^EPC|73000^X-RAY CLAVICLE^^^X-RAY CLAVIC|||7  
57^SPETZLER^CHRISTIAN^J^^^|(857)259-  
9167|||Final||^200310071730^^^|||2149^HITACHI^BENJAMIN^^^||1148010^1A  
^FRANKLIN^X-RAY^^^|^| InteleShare1||I9|810^CLAVICLE FRACTURE^I9|CLAVICLE  
FRACTURE||OBX|1|ST|GDT|1|X-ray suspected Clavicle fx. The quality of the film was  
good.|||A
```

```
|||Final|||1148010^1A^FRANKLIN^X-RAY^^^ OBX|2|ST|IMP|1| Abnormal, Fracture  
Present|||A|||Final|||1148010^1A^FRANKLIN^X-RAY^^^ OBX|3|ST|IMP|1|Left lateral;  
surgery required|||A|||Final|||1148010^1A^FRANKLIN^X-RAY^^^
```

ORM^O01 HL7 message

An ORM^O01 message is a general order message that contains ap- pointment and personal patient details. These are the incoming ORM workflows that InteleShare supports:

Outside study is sent to InteleShare

Outside orders are sent to InteleShare from the target facility. In InteleShare, the order is manually associated with study, and order information (MRN, accession number, patient name, DOB, gender, and study description) is used to update the study. After the update,

study is sent to PACS. Study flows automatically into PACS because it has the right target PHI.

ORMs are used to drive routing: sending to referring physicians, sharing with specified physicians, sending to DICOM devices, and any other DICOM destination.

The ORM auto associates with the study using accession number as a linking key, and the routing rule is configured to fire based on information in a specified segment and field.

Below is an example ORM^O01 message. The key field is the accession number. This is the linking key between the order and the study. The accession number is generally pulled from OBR1 or OBR2 (see red type in sample message below). The field must contain only the accession number. It cannot contain any ^ or sub-fields. Although OBR2 or OBR3 is typically used, any segment and field in the ORM^O01 message can be defined as field to link the report with the study.

See red text in sample message below.

```
MSH|^~\&|EMR|EMR|||20041111072701|1148|ORM^O01|497|D|2.3||
PID|1||43560^^^^EPI||ECKLES^SAMUEL^A.^MR.^||19730211|M||AfrAm|42 NORTH
AVE^^CAMBRI|InteleShareE^MA^02138^US^^DN|DN| DN|(857)721-3255-|(857)721-
3255|S||11480003~2514705~7603~13412|114-82-1544|||^^MA^^NTE|1||Comments
go here

PD1||FRANKLIN NEAR EAST LOCATION
(CENTRAL)^^537|757^SPETZLER^CHRISTIAN^J^^
PV1||^EMR^^^^|||509502|||V

ORC|NW|363463^EPC|1858^EPC||Final||^^200308151200^^|2004
01220727|1152^KARIS^LUCAS^^|757^SPETZLER^CHRISTIAN
^J^^|467^^537^^|(857)259-9167||

OBR|1|363463|1858^EPC|73610^X-RAY ANKLE 3+ VW^^X-RAY
ANKLE|||757^SPETZLER^CHRISTIAN^J^^(857)259-9167|||Final||^^200308
151200^^|2149^HITACHI^BENJAMIN^^|1148010^1A^FRANKLIN^X- RAY^^|^
InteleShare1||I9|824^ANKLE FRACTURE^I9|ANKLE FRACTURE||
```

ORMs are also used to drive DICOM modality worklist in the modality. The appointment/order is made in the EMR or practice management system; the order is sent to IntelShare, via IntelShare Gateway; IntelShare pulls the information from the ORM and acting as the service class provider (SCP) makes the information from the order available for DICOM modality worklist queries. When the patient arrives for a scan, instead of keying the MRN, accession number, name and other study information into the modality, the tech simply picks the patient name from a picklist.

IntelShare Gateway must be installed

IntelShare Gateway is the service class provider (SCP). This means that IntelShare Gateway receives MWL queries from modalities, which are acting as the service class user (SCU).

MWL is a standard supported by most modalities

For the modality to be set up to query the modality worklist, the customer should not need to pay the modality vendor. The customer should either be able to do it themselves or do it through basic modality support.

The queries are part of the log.log file in IntelShare Gateway

This is a useful resource during configuration.

Configuration

Appointment system/EMR/HIS

Configure the appointment system/EMR/HIS to send HL7 orders to IntelShare Gateway. You must provide the appointment system/EMR/HIS the following information:

- IP address of IntelShare Gateway
- HL7 listening port of IntelShare Gateway. By default, this is 12000. See hl7_listening_ports in IntelShare Gateway.
- The orders (ORM^O01) should contain the following information. Please note that fields formatted like OBR4-2 mean the second subfield in OBR4.
 - MRN in PID2 or PID3
 - Patient name in PID5

- Accession number in ORC3 or OBR3 or MHS10
- Birth date in PID7
- Sex in PID8
- Modality in OBR20
- Requested procedure id in OBR4-1 or ORC2
- Requested procedure description in OBR4-2
- Requested procedure step start date in OBR27-4
- Requested procedure station aetitle in OBR19

Below is an example ORM^O01 message that is part of the modality worklist workflow.

```
MSH|^~\&|EMR|EMR|||20041111072701|1148|ORM^O01|497|D|2.3||
PID|1||43560^^^^EPI||ECKLES^SAMUEL^A.^MR.^||19730211|M||AfrAm|4 2 NORTH
AVE^^CAMBRI|InteleShareE^MA^02138^US^^DN|DN| DN|(857)721-3255|(857)721-
3255||S||11480003~2514705~7603~13412|114-82-1544|||^^^MA^^
```

```
NTE|1||Comments go here
```

```
PD1|||FRANKLIN NEAR EAST LOCATION
(CENTRAL)^537|757^SPETZLER^CHRISTIAN^J^^PV1|||^EMR^^^^^|509
502|V
```

```
ORC|NW|363463^EPC|1858^EPC||Final||^200308151200^^^||200401220727|1
152^KARIS^LUCAS^^^||757^SPETZLER^CHRISTIAN^J^^|467^^^537^^^|
(857)259-9167||
```

```
OBR|1|363463|1858^EPC|73610^X-RAY ANKLE 3+ VW^^X-RAY ANKLE
```

```
|||||757^SPETZLER^CHRISTIAN^J^^|(857)259-9167|||||Final||^200308
151200^^^|||2149^HITACHI^BENJAMIN^^^||1148010^1A^FRANKLIN^X- RAY^^^|^|
```

```
InteleShare1||I9|824^ANKLE FRACTURE^I9|ANKLE FRACTURE||
```

InteleShare

Enable the orders_from_hl7 account switch. This must be done via the account/set call in the API.

Create the modality as a destination where:

- Port is filled in
- AE title is filled in
- IP address is filled in
- Type = DICOM (it defaults to DICOM if configured in the UI)

Remember that IntelShare is acting as the SCP in this process, meaning it is receiving queries from the modality.

Modality

Modality worklist needs to be configured on the modality. This typically means a configuration checkbox and these fields:

- SCP (IntelShare Gateway) IP address
- SCP (IntelShare Gateway) ae title
- SCP (IntelShare Gateway) port
- Also, please note that the modality must be able to talk with IntelShare Gateway.

Use

Patient arrives.

Tech picks patient from worklist on modality workstation. All information from appointment/order is written to study.

ADT^A08, ADT^A31, ADT^A04, ADT^A34 and SIU messages

ADT and SIU messages are used to sync data in your MPI (master patient index) and IntelShare. This is hugely valuable. This is how ADT and SIU messages are used in IntelShare:

Whenever there is a change to a patient record, an ADT message is sent to IntelShare.

The MRN is used to associate the ADT with the various studies that share the MRN. The PHI from the ADT including name, DOB, gender is written to the studies that share the MRN.

For ADT^A08 messages, MRN is updated through the value in PID2, name through the value in PID5, DOB through the value in PID7, and gender through the value in PID8.

For ADT^A34 messages, IntelShare looks at the change in value from MRG1 and updates using PID2, name through the value in PID5, DOB through the value in PID7, and gender through the value in PID8.

ADT Messages and SIU Messages can be used to create new patients.

ADT^A04 messages, ADT^A08, and SIU messages can be used to create new patients. For ADT^A04, ADT^A08, and SIU messages, Name almost always comes from PID5, DOB from PID7, and Gender from PID8. Email and phone depends on the version of HL7. Typically in pre-2.7, these segments are in PID13. In 2.7 or later, these segments are in PID40. IntelShare is capable of map- ping to whichever segments are needed for phone and email.

As a note, please leave PID19 (SSN) blank.

Below is an example ADT^A31 message. The key field is the MRN. This is the linking key between the ADT message and the study or studies. The MRN is generally pulled from PID2 or PID3.

See red text in sample message below.

```
MSH|^~\&|HQ|A^HH|||201506130206||ADT^A31|082057154L0027316745|D|2.2
EVN|A31|201508190206||HPIU|ANWILSON
PID|0001|90442421|090442421^^^HQA^MR~90442421^^^HQ^PI||SMITH^ALB
ERT||19920430|M||X||||ENG|S|CAT|||^01
PV1|0001|N|||||||||||||||||||||||||||||||||||||^PERSONALDOC
ZPI|0001|||||||20150818|N|||NONE^UNEMPLOYED
ZID|0001|||||04|""
ZV2||Y|||||||||U|||||N
ZEZ|N
```

7

OUTGOING HL7 MESSAGES

ORU^R01 HL7 messages

The IntelShare HL7 generator is very flexible. The messages described below are the most common messages IntelShare generates; however, IntelShare can generate other types of HL7 messages as required.

ORU^R01 messages are used to send reports to RIS, EMR or patient portal systems.

Common workflows include:

A report is made using IntelShare's reporting engine.

When complete, the report can be sent via ORU^R01 to a target RIS or EMR.

An encapsulated PDF is sent to an outside destination.

A PDF can be sent in an HL7 as a base64 encoded encapsulated PDF.

An encrypted URL link is commonly sent from IntelShare to an EMR or a patient portal via an ORU^R01 generated by IntelShare's HL7 engine.

For more information on encrypted URL links, please contact IntelShare. In addition to being an ORU generator, IntelShare can enrich ORUs received from other systems with an encrypted URL link or other data and then forward the enriched ORU to downstream RIS or EMR systems.

Below is an example ORU^R01 message with a URL link embedded in the message.

See red text in sample message below.

```
MSH|^~\&|RAD|^005^||||ORU|0000000741|P|2.5|
PID|||14055||SMITH^ANDREA||19520908|F||||||| PV1|||ENT|||||
OBR|||8999059211|NECK^1HEAD_NECK_PETCT|||||||
OBX|1|URL|||REDIRECT LINK GOES HERE|||||
```

ORM^O01 HL7 messages

IntelShare generates ORMs as part of our automated study PHI normalization process, which enables no touch outside study ingestion. The outside procedure code is automatically determined, the order is automatically created and written to the study, and the study is automatically pushed to PACS. The workflow is:

IntelShare receives an outside study that is being sent to a downstream PACS.

As the study is coming from the outside, the MRN, the accession number and possibly the patient name, date of birth and gender will need to be normalized to the target system's information. IntelShare generates an ORM and sends it to the target's RIS or EMR. The ORM includes the patient information, the outside procedure code, and all other relevant information that can be scraped off the study including modality, study description, study date, image count, and any other data the target RIS or EMR requires. The RIS or EMR receives IntelShare's order, automatically does a patient lookup, and if there is a match, automatically generates a real order, and sends the order to IntelShare. IntelShare receives the order, automatically applies the order information to the study, and approves the study, which triggers a send to PACS.

Below is an example ORM^O01 message that can help drive the workflow described above.

```
MSH|^~\&|SOARIAN|Target system|||20141119125329||ORM^O01|726976|T|2.4
PID|||1006551^^^MR ||ZZTRAIN^AMY||19630506000000|F|||||||
ORC|NW|
OBR||5789345|||||||O|FALSE|||
ZAI|DBID|||<Db ID>|||||||
```

8

MESSAGE SEGMENTS AND FIELDS

InteleShare supports message segments and fields as described here:
http://www.hl7.org/implement/standards/product_brief.cfm?product_id=185

9

RECENT UPDATES

Receiving HL7 Messages

InteleShare Gateway can be configured to receive HL7 messages in either UTF-8 (default) or Windows-1252 encoding.

HL7 Splitter

The HL7 splitter allows the Gateway to push HL7 messages directly to any HL7 destination on the local network. In parallel, the messages are uploaded to the InteleShare cloud. This does not require internet access, so if there is an outage, the gateway will continue to push messages to the local destination.

ADT Messages

You can now use ADT messages instead of ORMS to enable a modality worklist.

10

OBTAINING PRINTED DOCUMENTATION

Intelerad offers printed and bound versions of product documentation free of charge. To request printed copies of Intelerad documentation, contact your Client Success manager. The printed documents will be provided within 7 days or less.

11

CONTACTING INTELERAD TECHNICAL SUPPORT

Your PACS administrator can assist you with any issues you may encounter. If you require additional assistance, you can contact IntelShare Technical Support, 24 hours a day, seven days a week.

To contact us:	Use:
On the Internet	https://serviceportal.intelerad.com/csm
By telephone	Toll-free: 888-315-0790

These coordinates and a wealth of other information are also available on the Intelerad Service Portal.

<https://serviceportal.intelerad.com/csm>

You should regularly check the Intelerad knowledge base for the latest version of the documentation.

IntelShare PACS is a cloud-based product. Customers are informed of upcoming releases via multiple channels to ensure they stay up-to-date with the latest developments and enhancements. All customers receive timely and comprehensive e-mails from our marketing team, which not only highlight the release dates but also describe in detail the new features and functionalities that will be available in the upcoming release. When a new version of the software becomes available it will automatically upgrade on the appointed release date of that version. Please contact support with for any questions regarding upgrades.

When you contact IntelShare Technical Support to report a problem, please have at hand the following information, as applicable:

- client code and location of your IntelShare PACS installation
- full error message and the steps required to reproduce the problem
- AE Titles of the affected devices
- operating systems of any affected machines
- description of the problem and when it first occurred

If the problem affects a particular study, please also provide the following:

- patient ID or patient number (M.R.N.)
- accession number/requisition number
- modality type and name

For security information, please open a ticket with Intelrad Technical Support.